

**mRNA Cap 2'-O-Methyltransferase**

v. 240101

<b>Catalog number</b>	C15038-2000U	
<b>Package &amp; Component</b>	<b>Name</b>	<b>Amount</b>
	mRNA Cap 2'-O-Methyltransferase (50,000 U/mL)	2000 U (1 vial)
	10X Capping Enzyme Reaction Buffer	100 µL (1 vial)
<b>Description</b>	<p>mRNA Cap 2' O Methyltransferase specifically requires a Cap 0 structure (m7Gppp5'N) as a substrate, which synthesis from in vitro transcripts mRNA and capping enzyme modify.</p> <p>This enzyme utilizes a methyl donor such as SAM to add a methyl group at the 2' -O position forming cap-1 structure.</p>	
<b>Source</b>	<i>Escherichia coli</i>	
<b>Purity</b>	>95% as determined by SDS-PAGE. Purified by Ni-NTA chromatography.	
<b>Unit Definition</b>	One unit is defined as the amount of enzyme required to methylate 10 pmoles of 80 nt long capped RNA transcript in 1 hour at 37°C.	
<b>Reaction Condition</b>	<p>1X Capping enzyme reaction buffer, supplemented with 0.5 mM GTP and 0.1 mM S-adenosylmethionine (SAM). Incubate at 37°C.</p> <p>10X Capping enzyme Reaction Buffer: 500 mM Tris-HCl (pH 8.0), 50 mM KCl, 10 mM MgCl<sub>2</sub>, 10 mM DTT.</p>	
<b>Storage Buffer</b>	Capping enzyme is supplied in 20 mM Tris-HCl (pH 8.0), 100 mM NaCl, 1 mM DTT, 0.1 mM EDTA, 0.1% Triton® X-100 and 50% Glycerol.	
<b>Storage</b>	-20°C or -80°C for 12 months under sterile conditions from date of receipt.	

**One-Step Capping and 2'-O-Methylation procedures:**

1. Combine RNA and nuclease-free H<sub>2</sub>O to a final 14 µL.
2. Heating at 65°C for 5 minutes then chill on ice for 5 minutes.
3. Below reaction mixture should be prepared **on ice** and combined in the following order:

**Manuel**

<b>Component</b>	<b>Amount</b>	<b>Final concentration</b>
Denatured uncapped RNA	14 µL	-
10X Capping Buffer	2 µL	1X
10 mM GTP	1 µL	0.5 mM
4 mM SAM	1 µL	0.2 mM
Vaccinia Capping Enzyme	1 µL	10 U/rxn
mRNA Cap 2'-O-Methyltransferase	1 µL	50 U /rxn

4. Gently mix the reaction thoroughly to achieve uniform distribution.
5. Incubate at 37°C for 60 minutes (For RNA less than 200 nt long increase incubation time to 2 hours).

---

**Notes**

Reaction should be performed under RNase free condition. Use nuclease-free tubes, reagents and water to avoid RNase contamination. Also, wear gloves when working with RNA.

---

*For Research Use Only. Not for use in diagnostic and therapeutic procedures.*