


PRODUCT INFORMATION

Recombinant Der p 2 protein (mammalian cell)

v. 230301

Catalog number	C14045-100UG
Package	100 µg
Source	HEK293
Sequence	DQVDVKDCANHEIKKVLVPGCHGSEPCIIHRGKPFQLEAVFEANQNTKTAKIEIK ASIDGLEVDVPGIDPNACHYMKCPLVKGQQYDIKYTWNVPKIAPKSENVVVTVK VMGDDGVLACAIATHAKIRD with polyhistidine-SUMO tag at the N-terminus
Concentration	0.5 mg/mL
Purity	>90% as determined by SDS-PAGE.
Form	Liquid
Storage buffer	Phosphate buffered saline, pH 7.4.
Storage	Store at -20°C
Note	This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles.



SDS-PAGE analysis of Recombinant Der p 2 protein (mammalian cell)

For Research Use Only.