

PRODUCT INFORMATION

Monoclonal Anti-Bos d 7 Chimeric Human IgE, clone 1E11

v. 230301

Catalog number	C14008-100UG
Package	100 µg
Product type	Monoclonal, clone 1E11
Concentration	1 mg/mL
Immunogen	Bovine IgG
Isotype	Human ε, κ
Production	Hybridoma supernatant
Preparation	KappaSelect resin purified, 0.22 µm filtered
Storage buffer	Phosphate-buffered saline containing 0.09 % sodium azide, pH 7.2
Purity	> 95% by SDS-PAGE analysis
Storage	Store at -20°C. Avoid repeated freeze/thaw cycles. Suggestion: Aliquote the antibody into several vials.
Note	For research use only. Not for human in vivo or therapeutic use.

For Research Use Only.