

PRODUCT INFORMATION

Monoclonal Anti-Der p 1 Chimeric Human IgE, clone 20H4

v. 230301

Catalog number	C14013-100UG
Package	100 µg
Product type	Monoclonal, clone 20H4
Concentration	1 mg/mL
Immunogen	Recombinant Der p 1 fusion proteins
Specificity	Cross-react with recombinant Der f 1 fusion proteins
Isotype	Human ε, κ
Production	Mouse ascites
Preparation	KappaSelect resin purified, 0.22 µm filtered
Storage buffer	Phosphate-buffered saline containing 0.09 % sodium azide, pH 7.2
Purity	> 95% by SDS-PAGE analysis
Storage	Store at -20°C. Avoid repeated freeze/thaw cycles. Suggestion: Aliquote the antibody into several vials.
Note	Reactivity with natural or other recombinant Der p 1 or Der f 1 proteins is not tested. For research use only. Not for human in vivo or therapeutic use.

For Research Use Only.