

PRODUCT INFORMATION

Monoclonal Chimeric Human IgE Control, clone 29C11

v. 230301

Catalog number	C14028-100UG
Package	100 µg
Product type	Monoclonal, clone 29C11
Concentration	1 mg/mL
Immunogen	BSA-conjugated polypeptides
Specificity	Synthetic polypeptides
Isotype	Human ε, κ
Production	Mouse ascites
Preparation	KappaSelect resin purified, 0.22 µm filtered
Storage buffer	Phosphate-buffered saline containing 0.03 % Proclin 300, pH 7.4
Purity	> 95% by SDS-PAGE analysis
Storage	Store at -20°C. Avoid repeated freeze/thaw cycles. Suggestion: Aliquote the antibody into several vials.

For Research Use Only.