


## PRODUCT INFORMATION

**Recombinant Fel d 1 protein**

v. 230301

<b>Catalog number</b>	C14037-100UG
<b>Package</b>	100 µg
<b>Source</b>	<i>Escherichia coli</i>
<b>Sequence</b>	VKMAETCPIFYDVFFAVANGNELLLDLSLTKVNATEPERTAMKKIQDCYVENGLI SRVLDGLVMTTISSSKDCMGEAVQNTVEDLKLNTLGREICPAVKRDVDLFLTGT PDEYVEQVAQYKALPVVLENARILKNCVDAKMTEEDKENALSLLDKIYTSPLC with polyhistidine tag at the C-terminus
<b>Concentration</b>	0.5 mg/mL
<b>Purity</b>	>90% as determined by SDS-PAGE. Ni-NTA chromatography
<b>Form</b>	Liquid
<b>Storage buffer</b>	Phosphate buffered saline containing 0.01% sarkosyl, pH 7.4.
<b>Storage</b>	Store at -20°C
<b>Note</b>	This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles.



SDS-PAGE analysis of Recombinant Fel d 1 protein

*For Research Use Only.*