


PRODUCT INFORMATION

Recombinant Bos d 8 protein

v. 230301

Catalog number	C14040-100UG
Package	100 µg
Source	<i>Escherichia coli</i>
Sequence	RPKHPIKHQGLPQEVLNENLLRFFVAPFPEVFGKEKVNELSKDIGSESTEDQAM EDIKQMEAESISSSEEIVPNSVEQKHIQKEDVPSERYLGYLEQLLRLKKYKVPQL EIVPNSAEERLHSMKEGIHAQQKEPMIGVNQELAYFYPELFRQFYQLDAYPSGA WYYVPLGTQYTDAPSFSDIPNPIGSENSEKTTMPLW with polyhistidine tag at the C-terminus
Concentration	0.5 mg/mL
Purity	>85% as determined by SDS-PAGE. Ni-NTA chromatography
Form	Liquid
Storage buffer	Phosphate buffered saline containing 0.01% sarkosyl, pH 7.4.
Storage	Store at -20°C
Note	This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles.



SDS-PAGE analysis of Recombinant Bos d 8 protein

For Research Use Only.