

PRODUCT INFORMATION

Rabbit anti-Mouse IgG

v. 230201

Catalog number	C04044-1MG
Package	1 mg
Description	Croyez Rabbit anti-Mouse IgG is purified from goat antiserum by immunoaffinity chromatography. It reacts with mouse IgG heavy and light chains.
Product type	Polyclonal
Concentration	5 mg/mL
Immunogen	Purified mouse IgG
Host	Rabbit
Reactivity	Mouse
Conjugation	N/A
Isotype	N/A
Purity	Immunogen affinity chromatography.
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
Application	N/A
Application Note	N/A
Manual	N/A

For Research Use Only.