

PRODUCT INFORMATION

Mouse anti-Human IgG, clone hG15

v. 230201

Catalog number	C04047-1MG
Package	1 mg
Description	Croyez Mouse anti-Human IgG is a monoclonal antibody specific to human IgG.
Product type	Monoclonal, clone hG15
Concentration	1 mg/mL
Immunogen	Purified human IgG
Host	Mouse
Reactivity	Human
Conjugation	N/A
Isotype	IgG ₁ kappa
Purity	>95% (SDS-PAGE)
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 0.03% Proclin 300, pH 7.4.
Storage	Store at 4°C for at least two weeks. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Application	N/A
Application Note	N/A
Manual	N/A

For Research Use Only.