

## PRODUCT INFORMATION

**Mouse anti-Human IgG, clone hG17**

v. 230201

|                         |  |
|-------------------------|--|
| <b>Catalog number</b>   | C04048-1MG   |
| <b>Package</b>          | 1 mg   |
| <b>Description</b>      | Croyez Mouse anti-Human IgG is a monoclonal antibody specific to human IgG.                                |
| <b>Product type</b>     | Monoclonal, clone hG17   |
| <b>Concentration</b>    | 1 mg/mL  |
| <b>Immunogen</b>        | Purified human IgG   |
| <b>Host</b>             | Mouse  |
| <b>Reactivity</b>       | Human  |
| <b>Conjugation</b>      | N/A  |
| <b>Isotype</b>          | IgG <sub>1</sub> kappa   |
| <b>Purity</b>           | >95% (SDS-PAGE)  |
| <b>Form</b>             | Liquid   |
| <b>Storage buffer</b>   | Phosphate Buffered Saline containing 0.03% Proclin 300, pH 7.4.  |
| <b>Storage</b>          | Store at 4°C for at least two weeks.<br>Store at -20°C for one year.<br>Avoid repeated freeze/thaw cycles. |
| <b>Application</b>      | N/A  |
| <b>Application Note</b> | N/A  |
| <b>Manual</b>           | N/A  |

*For Research Use Only.*