

## PRODUCT INFORMATION

**Donkey anti-Goat IgG (H+L)-HRP**

v. 230201

<b>Catalog number</b>	C04006-1ML / C04006-2ML
<b>Package</b>	1 mL / 2 mL
<b>Description</b>	Croyez Donkey anti-Goat IgG (H+L)-HRP is conjugated with high-activity HRP. In direct ELISA, less than 1% cross-reactivity with mouse, rabbit and rat IgG is observed.
<b>Product type</b>	Polyclonal
<b>Immunogen</b>	Purified goat IgG
<b>Host</b>	Donkey
<b>Reactivity</b>	Goat (Minor cross-reactivity to human in higher concentration)
<b>Conjugation</b>	HRP (RZ>3)
<b>Isotype</b>	N/A
<b>Purity</b>	Immunogen affinity chromatography.
<b>Form</b>	Liquid
<b>Storage buffer</b>	Phosphate Buffered Saline containing 1% BSA, 0.01% thimerosal, and 50% glycerol, pH 7.4.
<b>Storage</b>	Store at 4°C for two weeks. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. Protect from light.
<b>Application</b>	ELISA, IHC, WB
<b>Application Note</b>	N/A
<b>Manual</b>	1:500-1:2500 for IHC 1:5000-1:25000 for ELISA and WB Working dilution for specific application should be determined by the investigator to obtain the best conditions.

*For Research Use Only.*