

PRODUCT INFORMATION

Donkey anti-Goat IgG

v. 230201

Catalog number	C04045-1MG
Package	1 mg
Description	Croyez Donkey anti-Goat IgG is purified from donkey antiserum by immunoaffinity chromatography. It reacts with goat IgG heavy and light chains.
Product type	Polyclonal
Concentration	Lot-dependent, please refer to product label.
Immunogen	Purified goat IgG
Host	Donkey
Reactivity	Goat (Minor cross-reactivity to human in higher concentration)
Conjugation	N/A
Isotype	N/A
Purity	Immunogen affinity chromatography.
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
Application	N/A
Application Note	N/A
Manual	N/A

For Research Use Only.