

PRODUCT INFORMATION

Normal Goat IgG

v. 230201

Catalog number	C06005-5MG / C06005-10MG
Package	5 mg / 10 mg
Description	Goat IgG is purified and isolated from pooled normal goat serum by protein A. The purified normal goat IgG may be used as an isotype control in various immunoassays, such as ELISA, WB, IP, FACS, ID and IE.
Product type	Polyclonal
Concentration	Lot-dependent, please refer to product label.
Immunogen	N/A
Host	Goat
Reactivity	N/A
Conjugation	N/A
Isotype	N/A
Purity	>95% (SDS-PAGE)
Form	Liquid
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
Application	ELISA, WB, IP, FACS, ID, IE
Application Note	N/A
Manual	N/A

For Research Use Only.