

## PRODUCT INFORMATION

**Mouse IgG3 Kappa Isotype Control Antibody**

v. 230201

<b>Catalog number</b>	C06009-bulk / C06009-1MG
<b>Package</b>	Customized package / 1 mg
<b>Description</b>	Mouse IgG3 Kappa isotype control antibody in used to estimate non-specific binding of target primary antibodies.
<b>Product type</b>	Monoclonal
<b>Concentration</b>	1 mg/mL
<b>Immunogen</b>	N/A
<b>Host</b>	Mouse
<b>Reactivity</b>	N/A
<b>Conjugation</b>	N/A
<b>Isotype</b>	IgG3 kappa
<b>Purity</b>	>95% (SDS-PAGE)
<b>Form</b>	Liquid
<b>Storage buffer</b>	Phosphate Buffered Saline containing 0.01% thimerosal, pH 7.4.
<b>Storage</b>	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles. Suggestion: Divide antibody into several vials. Keep only vials for usage at 4°C.
<b>Application</b>	ELISA, WB, IP, FACS, ID, IE
<b>Application Note</b>	N/A
<b>Manual</b>	N/A

*For Research Use Only.*