

PRODUCT INFORMATION

Protein G sepharose

v. 230201

Catalog number	C07006-250UL / C07006-2ML / C07006-5ML
Package	0.25 mL / 2 mL / 5 mL
Description	Recombinant protein G were conjugated to NHS-sepharose beads. Croyez Protein G sepharose has > 10 mg protein G conjugated per mL of settled sepharose.
Product type	Protein G
Concentration	N/A
Conjugation	Sepharose
Purity	N/A
Form	50% slurry
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C.
Application	IP
Application Note	N/A
Manual	The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of IgG protein. Optimal dilutions/concentrations should be determined by the end user.

For Research Use Only.