

PRODUCT INFORMATION

Ni-NTA sepharose

v. 230201

Catalog number	C07008-5ML / C07008-10ML / C07008-50ML
Package	5 mL / 10 mL / 50 mL
Description	Ni ²⁺ were conjugated to NHS-sepharose beads.
Concentration	N/A
Conjugation	Sepharose
Purity	N/A
Form	50% slurry
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C.
Application	Protein purification
Application Note	N/A
Manual	The binding and elution capacity of 1 mL sepharose are commonly more than 5 mg of recombinant 6X His fusion protein. Optimal dilutions/concentrations should be determined by the end user.

For Research Use Only.