

## PRODUCT INFORMATION

**Anti-EGFP sepharose**

v. 230201

<b>Catalog number</b>	C07009-250UL / C07009-2ML / C07009-5ML
<b>Package</b>	0.25 mL / 2 mL / 5 mL
<b>Description</b>	Anti-EGFP V <sub>H</sub> H antibody was conjugated to NHS-sepharose beads. Croyez Anti-EGFP sepharose has 6-10 mg mAb conjugated per mL of settled sepharose.
<b>Product type</b>	V <sub>H</sub> H
<b>Concentration</b>	N/A
<b>Immunogen</b>	Recombinant EGFP V <sub>H</sub> H
<b>Host</b>	N/A
<b>Conjugation</b>	Sepharose
<b>Isotype</b>	IgG2a kappa
<b>Purity</b>	N/A
<b>Form</b>	50% slurry
<b>Storage buffer</b>	20% ethanol
<b>Storage</b>	Store at 4°C.
<b>Application</b>	IP
<b>Application Note</b>	N/A
<b>Manual</b>	The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of GFP and EGFP fusion protein. Optimal dilutions/concentrations should be determined by the end user.

*For Research Use Only.*