

PRODUCT INFORMATION

Anti-mCherry sepharose

v. 230201

Catalog number	C07011-250UL / C07011-2ML / C07011-5ML
Package	0.25 mL / 2 mL / 5 mL
Description	Anti-mCherry V _H H antibody was conjugated to NHS-sepharose beads. Croyez Mouse anti-mcherry mAb sepharose has 6-10 mg mAb conjugated per mL of settled sepharose.
Product type	V _H H
Concentration	N/A
Immunogen	Recombinant mCherry protein
Host	N/A
Conjugation	Sepharose
Isotype	IgG2a kappa
Purity	N/A
Form	50% slurry
Storage buffer	20% ethanol
Storage	Store at 4°C.
Application	IP
Application Note	N/A
Manual	The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of RFP, mCherry or DsRed fusion protein. Optimal dilutions/concentrations should be determined by the end user.

For Research Use Only.