

PRODUCT INFORMATION

Mouse anti-DYKDDDDK mAb sepharose

v. 230201

Catalog number	C07014-250UL / C07014-2ML / C07014-5ML
Package	0.25 mL / 2 mL / 5 mL
Description	Mouse anti-DYKDDDDK monoclonal antibody were conjugated to NHS-sepharose beads. Croyez Mouse anti-DYKDDDDK mAb sepharose has 6-10 mg mAb conjugated per mL of settled sepharose.
Product type	Monoclonal
Concentration	N/A
Immunogen	Peptide: DYKDDDDK
Host	Mouse
Conjugation	Sepharose
Isotype	IgG2b kappa
Purity	N/A
Form	50% slurry
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage	Store at 4°C.
Application	IP
Application Note	N/A
Manual	The binding and elution capacity of 1 mL sepharose are commonly more than 1 mg of DYKDDDDK fusion protein. Optimal dilutions/concentrations should be determined by the end user.

For Research Use Only.