

## PRODUCT INFORMATION

**Mouse anti-human CD3 mAb**

v. 230601

<b>Catalog number</b>	C10002-bulk / C10002-100UG
<b>Package</b>	Customized package / 100 µg
<b>Description</b>	Mouse anti-human CD3 mAb recognize CD3 antigen of CD3/TCR complex found only on mature T cell or medullary thymocytes. This antibody is approved for the therapy of acute, glucocorticoid-resistant rejection of allogeneic renal, heart and liver transplants.
<b>Product type</b>	Monoclonal, clone OKT3
<b>Concentration</b>	1 mg/mL
<b>Immunogen</b>	N/A
<b>Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	N/A
<b>Isotype</b>	IgG2a
<b>Purity</b>	>95% (SDS-PAGE)
<b>Form</b>	Liquid
<b>Storage buffer</b>	Low endotoxin, azide free-phosphate buffered saline, 0.2 µm filter sterilized. Endotoxin level is less than 0.01 EU/µg of the antibody, as determined by the LAL test.
<b>Storage</b>	This product is stable after storage at: <ul style="list-style-type: none"><li>• 2-8°C for 6 months under sterile conditions from date of receipt.</li></ul> Briefly centrifuge vials before opening.
<b>Application</b>	IHC, FACS
<b>Application Note</b>	N/A

*For Research Use Only.*