

PRODUCT INFORMATION

Protein A sepharose

v. 240301

Catalog number	C07005-250UL / C07005-2ML / C07005-5ML
Package	0.25 mL / 2 mL / 5 mL
Description	Recombinant protein A were conjugated to NHS-sepharose beads. Croyez Protein A sepharose has > 10 mg protein A conjugated per mL of settled sepharose.
Product type	Protein A
Concentration	N/A
Conjugation	Sepharose
Purity	N/A
Feature	The binding capacity falls within the range of 15-20 mg of human or rabbit IgG/ mL of settled resin.Optimal dilutions/concentrations should be determined by the end user.
Form	50% slurry
Storage buffer	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
Storage & Stability	Store at 4°C for 1 year.
Application	IP
Application Note	N/A

For Research Use Only.